


| | | |
|---|---|--|
|  | Instruction | Réf. : IT-07 |
| | Medium (Z+Si+Vitamins) preparation | Version : 1 Date : 28/11/2011 Page : 1/3 |

Référence : based on Z-formula, of CHU-10 (UTCC), vitamines F/2 (UTCC)

Note : Use fresh Milli-Q water to prepare stock solutions and culture medium.
The stock solutions and culture medium are stored in the refrigerator.

List of stocks solutions :

- ZSV1 : Na_2SiO_3 , 5 H_2O at 21,2 g/L
- ZSV 2 : NaNO_3 at 85 g/L
- ZSV 3 : MgSO_4 , 7 H_2O at 25 g/L
- ZSV 4 : K_2HPO_4 , 3 H_2O at 40,62 g/L or (K_2HPO_4 at 20 g/L)
- ZSV 5 : $\text{Ca}(\text{NO}_3)_2$, 4 H_2O at 20 g/L
- ZSV 6 : Na_2CO_3 at 21 g/L
- ZSV 7 : Fe-EDTA solution (*)
- ZSV 8 : Trace elements solution (**)
- ZSV 9 : Vitamin B12 solution (*Cyanocobalamine*) at 0,0001 g/L (***)
- ZSV 10 : Biotin solution at 0,01 g/L
- ZSV 11 : Thiamine solution at 0,1 g/L

(*) Fe-EDTA solution:

(a) EDTA- Na_2 , 2 H_2O solution :

In a 500 ml volumetric flask, weigh 323,52 mg of *EDTA- Na_2 , 2 H_2O* and adjust to volume with 500 ml Milli-Q water.

(b) HCl 0,1M :


In a 100 ml volumetric flask, transfer 1 ml of concentrated HCl and adjust to volume with Milli-Q water.

(c) FeCl_3 , 6 H_2O solution :

In a 50 ml volumetric flask, weigh 1,351 g of FeCl_3 , 6 H_2O then adjust to volume with 0,1 M HCl.

In a 1 L volumetric flask, transfer 500 ml of solution (a) and 10 ml of solution (c)
Then adjust to volume with Milli-Q water.

| | Editor | Checker | Approving |
|------------|-----------------|---------------------|----------------|
| Name : | Huguet Isabelle | Chardon Cécile | Rimet Frédéric |
| Fonction : | Technician | Research technician | Ingineer |
| Visa : | | | |

| | | |
|---|---|--|
|  | Instruction | Réf. : IT-07 |
| | Medium (Z+Si+Vitamins) preparation | Version : 1 Date : 28/11/2011 Page : 2/3 |

(**) Trace elements solution :

In a 500 ml volumetric flask , weigh :

| Element name | Mass weight (en mg) |
|---|---------------------|
| MnSO ₄ , H ₂ O / MnSO ₄ , 4 H ₂ O | 858,5 / 1115 |
| Na ₂ WO ₄ , 2 H ₂ O | 16,5 |
| (NH ₄) ₆ Mo ₇ O ₂₄ , 4 H ₂ O | 44 |
| KBr | 59,5 |
| KI | 41,5 |
| ZnSO ₄ , 7 H ₂ O | 143,5 |
| Cd(NO ₃) ₂ , 4 H ₂ O | 77 |
| Co(NO ₃) ₂ , 6 H ₂ O | 73 |
| Cu SO ₄ , 5 H ₂ O | 62,5 |
| Ni(NO ₃) ₂ , 6 H ₂ O / NiSO ₄ (NH ₄) ₂ SO ₄ , 6 H ₂ O | 72,9 / 99 |
| Cr(NO ₃) ₃ , 9 H ₂ O / Cr(NO ₃) ₃ , 7 H ₂ O | 20,3 / 18,5 |
| NH ₄ VO ₃ / V ₂ O ₄ (SO ₄) ₃ , 16 H ₂ O | 5,5 / 17,5 |
| Al(SO ₄) ₂ K, 12 H ₂ O / Al ₂ (SO ₄) ₃ K ₂ SO ₄ , 24 H ₂ O | 237 / 237 |

Then adjust to volume with Milli-Q water.


(***) Vitamin B₁₂ solution (Cyanocobalamine) at 0,0001 g/L :

In a 100 ml volumetric flask, weigh 1 mg of vitamin B₁₂, then adjust to volume with Milli-Q water. This is the solution (a) .

Dilute 1ml of the solution (a) in a 100 ml volumetric flask with Milli-Q water.

This is the ZSV9 solution .

| | Editor | Checker | Approving |
|------------|-----------------|---------------------|----------------|
| Name : | Huguet Isabelle | Chardon Cécile | Rimet Frédéric |
| Fonction : | Technician | Research technician | Engineer |
| Visa : | | | |

| | | |
|---|---|--|
|  | Instruction | Réf. : IT-07 |
| | Medium (Z+Si+Vitamins) preparation | Version : 1 Date : 28/11/2011 Page : 3/3 |

Preparation of 1 liter of (Z+Si+Vitamins) culture medium:

| N_o of stock solution | Name of stock solution | Concentration of stock solution(g/L) | Volume of stock solution (mL) |
|--|---|---|--------------------------------------|
| ZSV1 | Na ₂ SiO ₃ , 5 H ₂ O | 21,2 | 2,65 |
| ZSV2 | NaNO ₃ | 85 | 5,50 |
| ZSV3 | MgSO ₄ , 7 H ₂ O | 25 | 1 |
| ZSV4 | K ₂ HPO ₄ , 3 H ₂ O | 40,62 | 1,2 |
| ZSV5 | Ca (NO ₃) ₂ , 4 H ₂ O | 20 | 2,95 |
| ZSV6 | Na ₂ CO ₃ | 21 | 1 |
| ZSV7 | Fe-EDTA solution | - | 10 |
| ZSV8 | Trace elements solution | - | 0,080 |
| ZSV9 | vitamin B12 solution | 0,0001 | 0,1 |
| ZSV10 | Biotin solution | 0,01 | 0,1 |
| ZSV11 | Thiamine solution | 0,1 | 2 |

Complete to 1 liter with with Milli-Q water.

Then filter the medium (Z+Si+Vitamins) through a filter of 0,22 µm in diameter in the laminar flow hood.

| | Editor | Checker | Approving |
|------------|-----------------|---------------------|----------------|
| Name : | Huguet Isabelle | Chardon Cécile | Rimet Frédéric |
| Fonction : | Technician | Research technician | Engineer |
| Visa : | | | |