

Note : This medium is used for the majority of chlorophyceae.  
 Use fresh Milli-Q water to prepare stock solutions and culture medium.  
 The stock solutions and culture medium are stored in the refrigerator.

**List of stocks solutions :**

- LC1 : Ca (NO<sub>3</sub>)<sub>2</sub>, 4 H<sub>2</sub>O at 20 g/L
- LC2 : K<sub>2</sub>HPO<sub>4</sub>, 3H<sub>2</sub>O at 40,62 g/L or (K<sub>2</sub>HPO<sub>4</sub> at 20 g/L)
- LC3 : MgSO<sub>4</sub>, 7 H<sub>2</sub>O at 25 g/L
- LC4 : NO<sub>3</sub>K at 50 g/L
- LC5 : Ferric ammonium citrate iron at 0,3 g/L
- LC6 : Soil extracts
- LC7 : Bryophytes extracts
- LC8 : Trace elements solution (\*)

(\*) Trace elements solution:

In a 500 ml volumetric flask, weigh :

Element name	Mass weight (en mg)
ZnCl <sub>2</sub>	52
MnCl <sub>2</sub> , 4 H <sub>2</sub> O	757
CuCl <sub>2</sub> , 2 H <sub>2</sub> O	26,5
CoCl <sub>2</sub> , 6 H <sub>2</sub> O	20
MoO <sub>4</sub> Na <sub>2</sub>	10,5

Then adjust to volume with Milli-Q water

	Editor	Checker	Approving
Name :	Huguet Isabelle	Chardon Cécile	Rimet Frédéric
Fonction :	Technician	Research technician	Engineer
Visa :			

**Preparation of 1 liter of L.C culture medium:**

<b>N<sub>o</sub> of stock solution</b>	<b>Name of stock solution</b>	<b>Concentration of stock solution(g/L)</b>	<b>Volume of stock solution (mL)</b>
LC1	Ca (NO <sub>3</sub> ) <sub>2</sub> , 4 H <sub>2</sub> O	20	1,5
LC2	K <sub>2</sub> HPO <sub>4</sub> , 3 H <sub>2</sub> O	40,62	1,30
LC3	MgSO <sub>4</sub> , 7 H <sub>2</sub> O	25	1,2
LC4	NO <sub>3</sub> K	50	4
LC5	Ferric ammonium citrate iron	0,3	10
LC6	Extract of Earth	-	10
LC7	Extract of foam	-	5
LC8	Trace elements solution	-	0,1

Complete to 1 liter with fresh MilliQ water.

Then filter the medium L-C through a filter of 0,22 µm in diameter in the laminar flow hood.

	<b>Editor</b>	<b>Checker</b>	<b>Approving</b>
Name :	Huguet Isabelle	Chardon Cécile	Rimet Frédéric
Fonction :	Technician	Research technician	Engineer
Visa :			